

Fi-E

Gutta Percha Obturation Device



Scan for more details in website

Wireless and electric

Perfect gutta-percha fluidity, uniform and smooth gutta-percha output without bubbles.

Easier than traditional manual gutta-percha injection to relieve hand fatigue for the operator.

Eliminate surgical hazards caused by gutta-percha leakage and overflow.

Maintenance-free thanks to automatic cleaning of remaining gutta-percha.

Reusable injection needle

360° rotatable

The gutta-percha bar is directly loaded into the injection needle.

There is no gutta-percha residue inside the machine, making it easier to clean.

The gutta-percha bar can be replaced at any time without waiting for the internal cooling of machine.



Dual button design

Symmetrically-designed dual injection buttons, sensitive and convenient.

Compact and lightweight Easier to hold

Diameter: 26mm*26mm

Length: 200mm

Weight: 148g

