Partnering with you to tackle your toughest material cutting challenges

Great White Bur 330

> Great White® Bur #2 (1558)



Restoration Removal - Cavity Preparation - Endodontic Access

simply better dentistry

irait Weite." er #1 (331)

SAVE 12 Hours per year, Reduce Overhead by 35%

Designed with technically advanced blade geometry and a thicker neck SS White[®] Great White[®] Burs are designed to eliminate bur breakage when cutting hard material such as metal, porcelain and composite.

"I have yet to experience anything that cuts as fast, accurate and smooth as a Great White" Bur, saving my practice time and money while reducing stress. A single Great White" Bur replaces my use of two to three carbide burs for PFM crown removal, a time savings of about 12 hours per year and over \$4,000 in productivity." - Howard S. Glazer, DDS

The Evolution of the GW1557 More Efficient Procedures by Design





Straight-bladed burs cut smoothly but are slower, especially in harder materials.

Dentate-bladed burs are faster cutting but have more potential for vibration.

Great White[®] **Burs** have a unique state-of -the-art-design. By increasing dentates and incorporating an enhanced blade geometry, the bur cuts more efficiently than a straight-bladed bur. They will not catch, grab or stall.



All-In-One Kit

Complete Restoration Removal Kit


Designed for the efficient atraumatic removal of failed amalgams, composites, and zirconia and PFM Crowns,

Kit Reorder No. 18180

- COMPLETE RESTORATION REMOVAL K

Kit Contents:

Great White[®] #1 Great White[®] #2 Revelation[®] Diamond 806-014C Revelation[®] Diamond 801-018C Revelation[®] Diamond 856-016C Great White[®] Z 801-018 Great White[®] Z 856-018



Multi-Functional Performance Reduces Overhead



Tested by Professionals...

"The Great White" #245 bur is an excellent instrument for generating conservative preparations, particularly in posterior teeth. It not only effectively cuts tooth structure but also amalgam and composite resin. Its geometric configuration is highly suitable for developing ideal cavity preparations for Class I and II composite resins."

- Karl Leinfelder DDS, MS

Endodontic Access



Crown Removal/Cutting



Cavity Preparations



Amalgam Removal



Frustrated with Breakage?

Product Strength (lbs)



A recent study at the University of Maryland confirms that SS White® 557 Burs are 84.67% stronger than Midwest® Standard 557 Burs.*

Not All Burs are Created Equal!

With SS White's industry-leading order fulfillment rate of 98%, and internal 9-point inspection process, dental dealers have stated that SS White® carbide burs have the lowest product return and customer complaint rate in the industry. **PRACTICE**



"SS White[®] carbide burs offer superior cutting ability, strength and consistency on all dental materials, eliminating the worry often associated with other carbide burs." - Louis Malcmacher DDS, MAGD

INSPIRATION