




# Free Scan Imaging Plate Scanner

**Guilin Woodpecker Medical Instrument Co., Ltd.**  
Information Industrial Park, Guilin National High-Tech  
Zone, Guilin, Guangxi, 541004 P.R. China  
Sales Dept: +86-773-5873196  
E-mail: woodpecker@glwoodpecker.com  
Website: <http://www.glwoodpecker.com>  
ZMN-XC-074



## Bringing you a fully-controlled and intelligent image viewing experience



### Specifications

Theoretical resolution: 17lp/mm  
Imaging plate sizes: 0, 1, 2, 3  
Connection mode: USB  
Grayscale: 65536 levels (16Bit)  
Working mode: Online mode

### Remarkable imaging sharpens each detail

25um scanning accuracy.  
17LP/mm+65535 grayscale detail recognition ability.  
Clear details and image processing make diagnosis and treatment more effective.




### Easy to operate



One-button start for scanning .  
120-second autolock  
to keep data safe.


### Advanced algorithm leads to intelligence

Multi-layered structural image algorithm framework  
enhances the details of super-resolution images.  
Multi-mode processing algorithm to meet the image viewing  
habits of different doctors.



### Quick image viewing

Imported ultra-thin imaging plate has  
moderate thickness to keep oral cavity  
comfortable.  
The 0-3 size imaging plates are available for  
photographing all tooth positions.





### Powerful chip Intelligent scanning

FPGA chip accelerates data acquisition.  
Support image viewing on mobile phones,  
tablets, and computers.  
Exclusive image viewing software.  
Greatly facilitate image viewing and processing.

### Real-time intelligent dose analysis Safe dental film photography

Appropriate dose.  
Accurately record the dose  
value for each photographing.  
Timely reminder of exposure  
dose for clear images.  
Ensure secure photographing.



### Good-looking and handy Compacted size of 0.046m<sup>2</sup> and space-saving position

The silver-black scheme is suitable  
for upscale and official occasions.  
The compacted size of 0.046 m<sup>2</sup> remains full-fledged  
Standing nearby armchair to safeguard your teeth

